



M18BP rabbit pAb

Catalog No	YP-Ab-09232
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MIS18BP1 C14orf106 KIAA1903 KNL2 M18BP1
Immunogen	Synthesized peptide derived from human M18BP AA range: 929-979
Specificity	This antibody detects endogenous levels of human M18BP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	M18BP;MIS18BP1 C14orf106 KIAA1903 KNL2 M18BP1;Mis18-binding protein 1 (Kinetochore-associated protein KNL-2 homolog) (HsKNL-2) (P243)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome, centromere . Associated with centromeres in interphase cells, from late anaphase to the G1 phase. Not detected on centromeres during earlier phases of mitosis. Associated with chromatin. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

