



LRC61 rabbit pAb

Catalog No	YP-Ab-09233
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LRRC61
Immunogen	Synthesized peptide derived from human LRC61 AA range: 94-144
Specificity	This antibody detects endogenous levels of human LRC61
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRC61;LRRC61;Leucine-rich repeat-containing protein 61
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

