



K132L rabbit pAb

Catalog No	YP-Ab-09235
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KIAA1324L EIG121L
Immunogen	Synthesized peptide derived from human K132L AA range: 735-785
Specificity	This antibody detects endogenous levels of human K132L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	K132L;KIAA1324L EIG121L;UPF0577 protein KIAA1324-like (Estrogen-induced gene 121-like protein) (hEIG121L)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

