



CLCA rabbit pAb

Catalog No	YP-Ab-09237
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	CLTA
Immunogen	Synthesized peptide derived from human CLCA AA range: 188-238
Specificity	This antibody detects endogenous levels of human CLCA
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CaCC;CaCC-1;CACCC1;Calcium dependent chloride channel 1;Calcium-activated chloride channel family member 1;Calcium-activated chloride channel protein 1;Calcium-activated chloride channel regulator 1;Chloride channel calcium activated family member 1;Chloride channel regulator 1;CLCA family member 1 chloride channel regulator;CLCA1;CLCA1_HUMAN;CLCRG1;GOB5;hCaCC-1;hCLCA1;CLCA;CLTA;Clathrin light chain A (Lca)
Calculated Molecular Weight	1kDa; 1kD
Cell Pathway	Cytoplasmic vesicle membrane; Peripheral membrane protein; Cytoplasmic side. Membrane, coated pit; Peripheral membrane protein; Cytoplasmic side. Cytoplasm, cytoskeleton, spindle . Cytoplasmic face of coated pits and vesicles. In complex with TACC3 and CKAP5 (forming the TACC3/ch-TOG/clathrin complex) localized to inter-microtubule bridges in mitotic spindles. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images