



NOE1 rabbit pAb

Catalog No	YP-Ab-09239
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	OLFM1 NOE1 NOEL1
Immunogen	Synthesized peptide derived from human NOE1 AA range: 414-464
Specificity	This antibody detects endogenous levels of human NOE1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AMY;NEUROBLASTOMA PROTEIN;Neuronal olfactomedin-related ER localized protein;Noe1;NOE1_HUMAN;Noel;Noel1;Noelin;Noelin1;OlfA; Olfactomedin 1;Olfactomedin-1;Olfm1;Pancortin 1;Pancortin;OLFM1 NOE1 NOEL1;Noelin (Neuronal olfactomedin-related ER localized protein) (Olfactomedin-1)
Observed Band	
Calculated Molecular Weight	55kDa; 55kD
Cell Pathway	Secreted . Cell junction, synapse . Endoplasmic reticulum . Cell projection, axon . Perikaryon .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

