



TM108 rabbit pAb

Catalog No	YP-Ab-09240
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TMEM108 KIAA1690 UNQ1875/PRO4318
Immunogen	Synthesized peptide derived from human TM108 AA range: 88-138
Specificity	This antibody detects endogenous levels of human TM108
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM108;TMEM108 KIAA1690 UNQ1875/PRO4318;Transmembrane protein 108
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein . Cell junction, synapse, postsynaptic density . Endosome membrane . Cell projection, axon . Cell projection, dendrite . Early endosome .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

