



LYZL2 rabbit pAb

Catalog No	YP-Ab-09241
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	LYZL2
Immunogen	Synthesized peptide derived from human LYZL2 AA range: 22-72
Specificity	This antibody detects endogenous levels of human LYZL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LYZL2;Lysozyme-like protein 2 (Lysozyme-2) (EC 3.2.1.17)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in testis, epididymis and placenta.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

