



E4F1 rabbit pAb

Catalog No	YP-Ab-09242
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	E4F1 E4F
Immunogen	Synthesized peptide derived from human E4F1
Specificity	This antibody detects endogenous levels of human E4F1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E4F;E4F transcription factor 1;isoform CRA_b;p12E4F;P12E4F transcription factor;p5E4F;Putative E3 ubiquitin-protein ligase E4F1;RING type E3 ubiquitin transferase E4F1;Transcription factor E4F1;E4F1;E4F1 E4F;Transcription factor E4F1 (EC 6.3.2.-) (E4F transcription factor 1) (Putative E3 ubiquitin-protein ligase E4F1) (Transcription factor E4F) (p12E4F) (p50E4F)
Observed Band	
Calculated Molecular Weight	83kD
Cell Pathway	Nucleus, nucleoplasm. Cytoplasm. A small fraction is detected in the cytoplasm. Excluded from the nucleolus where it is targeted upon CDKN2A overexpression. Localizes to the mitotic spindle during embryogenesis. .
Tissue Specificity	Ubiquitously expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

