



YN4403

LHX6 rabbit pAb

Catalog No	YP-Ab-09243
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LHX6 LHX6.1
Immunogen	Synthesized peptide derived from human LHX6 AA range: 257-307
Specificity	This antibody detects endogenous levels of human LHX6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LHX6.1;LIM homeobox 6;LIM homeobox protein 6;LIM/homeobox protein Lhx6; LIM/homeobox protein Lhx6.1;LHX6;LHX6 LHX6.1;LIM/homeobox protein Lhx6 (LIM homeobox protein 6) (LIM/homeobox protein Lhx6.1)
Observed Band	40 kDa
Calculated Molecular Weight	392 aa, 43 kDa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images