



HN1L rabbit pAb

Catalog No	YP-Ab-09244
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HN1L C16orf34 L11
Immunogen	Synthesized peptide derived from human HN1L AA range: 135-185
Specificity	This antibody detects endogenous levels of human HN1L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C16orf34;Hematological and neurological expressed 1-like protein;HN1-like protein;JPT2;Jupiter microtubule associated homolog 2;HN1L;HN1L C16orf34 L11;Hematological and neurological expressed 1-like protein (HN1-like protein)
Observed Band	20-25 kDa
Calculated Molecular Weight	190 aa, 20 kDa
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Expressed in liver, kidney, prostate, testis and uterus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

