



CHM2B rabbit pAb

Catalog No	YP-Ab-09248
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	CHMP2B CGI-84
Immunogen	Synthesized peptide derived from human CHM2B AA range: 53-103
Specificity	This antibody detects endogenous levels of human CHM2B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALS17;Charged multivesicular body protein 2b;CHM2B_HUMAN;CHMP family; member 2B;CHMP2.5;CHMP2b;Chromatin modifying protein 2b;Chromatin-modifying protein 2b;DMT1;hVps2-2;Vacuolar protein sorting 2;yeast;homolog of; B;Vacuolar protein sorting 2-2;Vacuolar protein sorting-associated protein 2-2; VPS2 homolog B;Vps2-2;VPS2B;CHM2B;CHMP2B CGI-84;Charged multivesicular body protein 2b (CHMP2.5) (Chromatin-modifying protein 2b) (CHMP2b) (Vacuolar protein sorting-associated protein 2-2) (Vps2-2) (hVps2-2)
Observed Band	
Calculated Molecular Weight	32 kDa; 24kD
Cell Pathway	Cytoplasm, cytosol . Late endosome membrane ; Peripheral membrane protein .
Tissue Specificity	Widely expressed. Expressed in brain, heart, skeletal muscle, spleen, kidney, liver, small intestine, pancreas, lung, placenta and leukocytes. In brain, it is expressed in cerebellum, cerebral cortex, medulla, spinal chord, occipital lobe, frontal lobe, temporal lobe and putamen.
Function	
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images