



RF1ML rabbit pAb

Catalog No	YP-Ab-09249
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	MTRF1L MTRF1A
Immunogen	Synthesized peptide derived from human RF1ML AA range: 301-351
Specificity	This antibody detects endogenous levels of human RF1ML
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RF1ML;MTRF1L MTRF1A;Peptide chain release factor 1-like;mitochondrial (Mitochondrial translational release factor 1-like) (mtRF1a)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion .
Tissue Specificity	Expressed in skeletal muscle (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

