



DHX35 rabbit pAb

Catalog No	YP-Ab-09250
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DHX35 C20orf15 DDX35
Immunogen	Synthesized peptide derived from human DHX35 AA range: 618-668
Specificity	This antibody detects endogenous levels of human DHX35
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C20orf15;DDX35;DEAH box protein 35;EC:3.6.4.13;KAIA0875;DHX35;DHX35 C20orf15 DDX35;Probable ATP-dependent RNA helicase DHX35 (EC 3.6.4.13) (DEAH box protein 35)
Observed Band	70 kDa
Calculated Molecular Weight	703 aa, 79 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

