



YIF1A rabbit pAb

Catalog No	YP-Ab-09254
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	YIF1A 54TM HYIF1P YIF1
Immunogen	Synthesized peptide derived from human YIF1A AA range: 101-151
Specificity	This antibody detects endogenous levels of human YIF1A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	YIF1A;YIF1A 54TM HYIF1P YIF1;Protein YIF1A (54TMp) (YIP1-interacting factor homolog A)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Golgi apparatus membrane ; Multi-pass membrane protein . Endoplasmic reticulum-Golgi intermediate compartment membrane ; Multi-pass membrane protein . Cycles between the endoplasmic reticulum and the endoplasmic reticulum-Golgi intermediate compartment. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images