



MB3L1 rabbit pAb

Catalog No	YP-Ab-09257
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	MBD3L1 MBD3L
Immunogen	Synthesized peptide derived from human MB3L1 AA range: 136-186
Specificity	This antibody detects endogenous levels of human MB3L1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MB3L1;MBD3L1 MBD3L;Methyl-CpG-binding domain protein 3-like 1 (MBD3-like protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nuclear, in large foci.
Tissue Specificity	Highly expressed in testis. Detected at low levels in pancreas. Not detected in the other tissues tested.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

