



SAMD9 rabbit pAb

Catalog No	YP-Ab-09258
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SAMD9 C7orf5 DRIF1 KIAA2004 OEF1
Immunogen	Synthesized peptide derived from human SAMD9 AA range: 499-549
Specificity	This antibody detects endogenous levels of human SAMD9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C7orf5;DRIF1;KIAA2004;NFTC;OEF1;SAMD9;SAMD9 C7orf5 DRIF1 KIAA2004 OEF1;Sterile alpha motif domain-containing protein 9 (SAM domain-containing protein 9)
Observed Band	148 kDa
Calculated Molecular Weight	1589 aa, 184 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Widely expressed. Very low levels are detected in skeletal muscle. Not detected in brain. Down-regulated in aggressive fibromatosis, as well as in breast and colon cancers. Up-regulated in fibroblasts from patients with normophosphatemic tumoral calcinosis (NFTC).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images