



SYNG1 rabbit pAb

Catalog No	YP-Ab-09260
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	SYNDIG1 C20orf39 TMEM90B
Immunogen	Synthesized peptide derived from human SYNG1 AA range: 195-245
Specificity	This antibody detects endogenous levels of human SYNG1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SYNG1;SYNDIG1 C20orf39 TMEM90B;Synapse differentiation-inducing gene protein 1 (SynDIG1) (Dispanin subfamily C member 2) (DSPC2) (Transmembrane protein 90B)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane; Single-pass type II membrane protein. Early endosome membrane ; Single-pass type II membrane protein . Cell junction, synapse, postsynaptic density membrane . Cell junction, synapse . Cell projection, dendrite . Cell projection, dendritic spine . Shuttles between the cell surface and early endosome membrane. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images