



REBL1 rabbit pAb

Catalog No	YP-Ab-09262
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RHEBL1
Immunogen	Synthesized peptide derived from human REBL1 AA range: 32-82
Specificity	This antibody detects endogenous levels of human REBL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Activating transcription factor 6 beta;ATF 6 beta;ATF6-beta;ATF6B; ATF6B_HUMAN;ATF6beta;cAMP response element binding protein related protein;cAMP response element-binding protein-related protein;cAMP responsive element binding protein like 1;cAMP-dependent transcription factor ATF-6 beta;cAMP-responsive element-binding protein-like 1;CREB L1;Creb related protein;CREB RP;Creb-rp;CREBL 1;CREBL1;Crebrp;Cyclic AMP dependent transcription factor ATF 6 beta;Cyclic AMP dependent transcription factor ATF6 beta;FLJ166;G13;G13 protein;OTTHUMP29432;OTTHUMP29433; Processed cyclic AMP-dependent transcription factor ATF-6 beta;Protein G13; REBL1;RHEBL1;GTPase RhebL1 (Ras homolog enriched in brain like-1 c) (RhebL1c) (Ras homolog enriched in brain-like protein 1) (Rheb-like protein 1) (Rheb2)
Observed Band	
Calculated Molecular Weight	76 KD; 77kD
Cell Pathway	Endomembrane system ; Lipid-anchor ; Cytoplasmic side . Cytoplasm .
Tissue Specificity	Ubiquitously expressed. Expression increased at least 2-fold in several tumor cell lines.
Function	



Background

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

