



TM10B rabbit pAb

Catalog No	YP-Ab-09263
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TRMT10B RG9MTD3
Immunogen	Synthesized peptide derived from human TM10B AA range: 120-170
Specificity	This antibody detects endogenous levels of human TM10B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM10B;TRMT10B RG9MTD3;tRNA methyltransferase 10 homolog B (EC 2.1.1.-) (RNA (guanine-9-)-methyltransferase domain-containing protein 3)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

