



DPH3B rabbit pAb

Catalog No	YP-Ab-09265
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	DPH3P1 C20orf143 DPH3B ZCSL1
Immunogen	Synthesized peptide derived from human DPH3B AA range: 23-73
Specificity	This antibody detects endogenous levels of human DPH3B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DPH3B;DPH3P1 C20orf143 DPH3B ZCSL1;Putative DPH3 homolog B (CSL-type zinc finger-containing protein 1)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

