



UTY rabbit pAb

Catalog No	YP-Ab-09267
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	UTY
Immunogen	Synthesized peptide derived from human UTY AA range: 1096-1146
Specificity	This antibody detects endogenous levels of human UTY
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Histone demethylase UTY;KDM6AL;KDM6C;UTY1;UTY;Histone demethylase UTY (EC 1.14.11.-) (Ubiquitously-transcribed TPR protein on the Y chromosome) (Ubiquitously-transcribed Y chromosome tetratricopeptide repeat protein)
Observed Band	130-160kDa
Calculated Molecular Weight	150kDa,1347aa
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

