



YN4429

K1468 rabbit pAb

Catalog No	YP-Ab-09269
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	KIAA1468
Immunogen	Synthesized peptide derived from human K1468 AA range: 525-575
Specificity	This antibody detects endogenous levels of human K1468
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	K1468;KIAA1468;LisH domain and HEAT repeat-containing protein KIAA1468
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Recycling endosome . Golgi apparatus, trans-Golgi network . Translocated to the trans-Golgi network area in an OSBP-dependent manner (PubMed:29514919). Colocalizes with RAB11A in recycling endosomes (By similarity). Found in a complex composed of RELCH, OSBP1 and RAB11A (By similarity). .
Tissue Specificity	
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images