



CL029 rabbit pAb

Catalog No	YP-Ab-09270
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	C12orf29
Immunogen	Synthesized peptide derived from human CL029 AA range: 254-304
Specificity	This antibody detects endogenous levels of human CL029
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C12orf29;Chromosome 12 open reading frame 29;CL29_HUMAN;LOC9128;MGC12978;Uncharacterized protein C12orf29;CL029
Observed Band	
Calculated Molecular Weight	4,45kD; 37kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

