



LYZL4 rabbit pAb

Catalog No	YP-Ab-09271
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LYZL4 LYC4
Immunogen	Synthesized peptide derived from human LYZL4 AA range: 88-138
Specificity	This antibody detects endogenous levels of human LYZL4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	1810009N24Rik;LYC4;Lysozyme 4;lysozyme like 4;Lysozyme like protein 4; LYZL4;LYZL4 LYC4;Lysozyme-like protein 4 (Lysozyme-4)
Observed Band	16 kDa
Calculated Molecular Weight	146 aa, 16 kDa
Cell Pathway	Secreted . Cytoplasmic vesicle, secretory vesicle, acrosome . Cell projection, cilium, flagellum . Found in the principal piece of sperm tail. .
Tissue Specificity	Expressed in testis and epididymis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

