



STPG2 rabbit pAb

Catalog No	YP-Ab-09272
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	STPG2 C4orf37
Immunogen	Synthesized peptide derived from human STPG2 AA range: 102-152
Specificity	This antibody detects endogenous levels of human STPG2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	STPG2;STPG2 C4orf37;Sperm-tail PG-rich repeat-containing protein 2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

