



TM110 rabbit pAb

Catalog No	YP-Ab-09273
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	TMEM110
Immunogen	Synthesized peptide derived from human TM110 AA range: 220-270
Specificity	This antibody detects endogenous levels of human TM110
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM110;TMEM110;Transmembrane protein 110
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein . Colocalizes with STIM1 at ER-plasma membrane (ER-PM) junctions, also called cortical endoplasmic reticulum (ER), in store-depleted calcium cells (PubMed:26644574) . May translocate to ER-PM junctions in a STIM1-dependent manner in store-depleted cells (Probable). .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images