



ARF3 rabbit pAb

Catalog No	YP-Ab-09276
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ARF3
Immunogen	Synthesized peptide derived from human ARF3 AA range: 66-116
Specificity	This antibody detects endogenous levels of human ARF3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ADP ribosylation factor 3;ADP-ribosylation factor 3;ARF3
Observed Band	18-20 kDa
Calculated Molecular Weight	21 kDa
Cell Pathway	Golgi apparatus . Cytoplasm, perinuclear region .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

