



RL23A rabbit pAb

Catalog No	YP-Ab-09278
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	RPL23A
Immunogen	Synthesized peptide derived from human RL23A AA range: 17-67
Specificity	This antibody detects endogenous levels of human RL23A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	6S ribosomal protein L23A;FLJ27455;L23A;MDA2;Melanoma differentiation-associated gene 2;Ribosomal protein L23a;RL23A_HUMAN;RPL23A;RL23A;60S ribosomal protein L23a
Observed Band	
Calculated Molecular Weight	18kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

