



TTC31 rabbit pAb

Catalog No	YP-Ab-09280
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TTC31
Immunogen	Synthesized peptide derived from human TTC31 AA range: 419-469
Specificity	This antibody detects endogenous levels of human TTC31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TPR repeat protein 31;TTC31;Tetratricopeptide repeat protein 31 (TPR repeat protein 31)
Observed Band	57 kDa
Calculated Molecular Weight	358 aa, 40 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

