



SEHL2 rabbit pAb

Catalog No	YP-Ab-09282
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SERHL2 SERHL
Immunogen	Synthesized peptide derived from human SEHL2 AA range: 259-309
Specificity	This antibody detects endogenous levels of human SEHL2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SEHL2;SERHL2 SERHL;Serine hydrolase-like protein 2 (EC 3.1.-.-)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, perinuclear region . Peroxisome . Concentrated in perinuclear vesicles. May be located in peroxisomes. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

