



HMN9 rabbit pAb

Catalog No	YP-Ab-09284
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	MTRNR2L9
Immunogen	Synthesized peptide derived from human HMN9 AA range: 1-24
Specificity	This antibody detects endogenous levels of human HMN9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HMN9;MTRNR2L9;Humanin-like protein 9 (HN9) (MT-RNR2-like protein 9)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted . Cytoplasm .
Tissue Specificity	Highly expressed in the kidney, heart muscle and testis.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

