



LAD1 rabbit pAb

Catalog No	YP-Ab-09286
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LAD1 LAD
Immunogen	Synthesized peptide derived from human LAD1 AA range: 404-454
Specificity	This antibody detects endogenous levels of human LAD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LAD;Lad 1;LAD1;LadA;ladinin 1;Linear IgA disease antigen;LAD1 LAD;Ladinin-1 (Lad-1) (Linear IgA disease antigen) (LADA)
Observed Band	57 kDa
Calculated Molecular Weight	517 aa, 57 kDa
Cell Pathway	Secreted, extracellular space, extracellular matrix, basement membrane .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

