



ST65G rabbit pAb

Catalog No	YP-Ab-09287
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SUPT7L KIAA0764
Immunogen	Synthesized peptide derived from human ST65G AA range: 191-241
Specificity	This antibody detects endogenous levels of human ST65G
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ST65G;SUPT7L KIAA0764;STAGA complex 65 subunit gamma (Adenocarcinoma antigen ART1) (SPTF-associated factor 65 gamma) (STAF65gamma) (Suppressor of Ty 7-like)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed at high levels in adenocarcinomas and gliomas and low in esophageal cancers and malignant hematological disease. Also expressed at high level in the thymus, low in peripheral blood mononuclear cells, and lowest in the stomach, small intestine, and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images