



CEL3A rabbit pAb

Catalog No	YP-Ab-09288
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CELA3A ELA3 ELA3A
Immunogen	Synthesized peptide derived from human CEL3A AA range: 29-79
Specificity	This antibody detects endogenous levels of human CEL3A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CEL3A;CELA3A ELA3 ELA3A;Chymotrypsin-like elastase family member 3A (EC 3.4.21.70) (Elastase IIIA) (Elastase-3A) (Protease E)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

