



ATMIN rabbit pAb

Catalog No	YP-Ab-09289
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	ATMIN KIAA0431 ZNF822
Immunogen	Synthesized peptide derived from human ATMIN AA range: 378-428
Specificity	This antibody detects endogenous levels of human ATMIN
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATMIN;ATMIN KIAA0431 ZNF822;ATM interactor (ATM/ATR-substrate CHK2-interacting zinc finger protein) (ASCIZ) (Zinc finger protein 822)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Nuclear, in discrete foci during G1 phase.
Tissue Specificity	Ubiquitously expressed in normal tissues and cancer cell lines with highest levels in placenta and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

