



TECT1 rabbit pAb

Catalog No	YP-Ab-09290
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	TCTN1 TECT1 UNQ9369/PRO34160
Immunogen	Synthesized peptide derived from human TECT1 AA range: 68-118
Specificity	This antibody detects endogenous levels of human TECT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UNQ9369/PRO34160;Tectonic-1;Tectonic 1;TECT1;TCTN1 TECT1 UNQ9369/PRO34160
Observed Band	55-64 kDa
Calculated Molecular Weight	64 kDa
Cell Pathway	Cytoplasm, cytoskeleton, cilium basal body . Secreted . Despite the presence of a signal sequence, the full-length protein might not be secreted. Localizes at the transition zone, a region between the basal body and the ciliary axoneme. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

