



PGRP2 rabbit pAb

Catalog No	YP-Ab-09292
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	PGLYRP2 PGLYRPL PGRPL UNQ3103/PRO10102
Immunogen	Synthesized peptide derived from human PGRP2 AA range: 90-140
Specificity	This antibody detects endogenous levels of human PGRP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PGRP2;PGLYRP2 PGLYRPL PGRPL UNQ3103/PRO10102;N-acetylmuramoyl-L-alanine amidase (EC 3.5.1.28) (Peptidoglycan recognition protein 2) (Peptidoglycan recognition protein long) (PGRP-L)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted. Membrane.
Tissue Specificity	Strongly expressed in liver and fetal liver, and secreted into serum. Expressed to a much lesser extent in transverse colon, lymph nodes, heart, thymus, pancreas, descending colon, stomach and testis. Isoform 2 is not detected in the liver or serum.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images