



ZN483 rabbit pAb

Catalog No	YP-Ab-09294
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF483 KIAA1962 ZKSCAN16
Immunogen	Synthesized peptide derived from human ZN483 AA range: 315-365
Specificity	This antibody detects endogenous levels of human ZN483
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN483;ZNF483 KIAA1962 ZKSCAN16;Zinc finger protein 483 (Zinc finger protein with KRAB and SCAN domains 16)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

