



RD21L rabbit pAb

Catalog No	YP-Ab-09297
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RAD21L1 RAD21L
Immunogen	Synthesized peptide derived from human RD21L AA range: 501-551
Specificity	This antibody detects endogenous levels of human RD21L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RD21L;RAD21L1 RAD21L;Double-strand-break repair protein rad21-like protein 1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Chromosome . In meiotic chromosomes, localized along axial elements in early meiosis: detectable on the axial elements in leptotene, and stays on the axial/lateral elements until mid pachytene. It then disappears and is replaced with RAD21. Compared to REC8, has mutually exclusive loading sites on the chromosomes: REC8 and RAD21L form distinct cohesin-enriched domains along the axial elements. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images