



LRCH4 rabbit pAb

Catalog No	YP-Ab-09300
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LRCH4 LRN LRRN1 LRRN4
Immunogen	Synthesized peptide derived from human LRCH4 AA range: 145-195
Specificity	This antibody detects endogenous levels of human LRCH4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRCH4;LRCH4 LRN LRRN1 LRRN4;Leucine-rich repeat and calponin homology domain-containing protein 4 (Leucine-rich repeat neuronal protein 4) (Leucine-rich neuronal protein)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

