



H2A2C rabbit pAb

Catalog No	YP-Ab-09301
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HIST2H2AC H2AFQ
Immunogen	Synthesized peptide derived from human H2A2C AA range: 70-120
Specificity	This antibody detects endogenous levels of human H2A2C
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	H2A2C;HIST2H2AC H2AFQ;Histone H2A type 2-C (Histone H2A-GL101) (Histone H2A/q)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus. Chromosome.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

