



U3IP2 rabbit pAb

Catalog No	YP-Ab-09302
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	RRP9 RNU3IP2 U355K
Immunogen	Synthesized peptide derived from human U3IP2 AA range: 198-248
Specificity	This antibody detects endogenous levels of human U3IP2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	U3IP2;RRP9 RNU3IP2 U355K;U3 small nucleolar RNA-interacting protein 2 (RRP9 homolog) (U3 small nucleolar ribonucleoprotein-associated 55 kDa protein) (U3 snoRNP-associated 55 kDa protein) (U3-55K)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

