



UK114 rabbit pAb

Catalog No	YP-Ab-09303
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HRSP12 PSP
Immunogen	Synthesized peptide derived from human UK114 AA range: 52-102
Specificity	This antibody detects endogenous levels of human UK114
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	heat responsive protein 12;HRSP12;P14.5;PSP;Ribonuclease UK114;UK114; UK114 antigen homolog;HRSP12 PSP;Ribonuclease UK114 (EC 3.1.-.-) (14.5 kDa translational inhibitor protein) (p14.5) (Heat-responsive protein 12) (UK114 antigen homolog)
Observed Band	14.5 kDa
Calculated Molecular Weight	14.5 kDa
Cell Pathway	Cytoplasm . Nucleus . Peroxisome . Mitochondrion . Mostly cytoplasmic but, in less differentiated cells occasionally nuclear. .
Tissue Specificity	Expressed predominantly in liver and kidney. Lower levels in lung and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

