



S23IP rabbit pAb

Catalog No	YP-Ab-09304
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	SEC23IP MSTP053
Immunogen	Synthesized peptide derived from human S23IP AA range: 823-873
Specificity	This antibody detects endogenous levels of human S23IP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	S23IP;SEC23IP MSTP053;SEC23-interacting protein (p125)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle, COPII-coated vesicle membrane; Peripheral membrane protein; Cytoplasmic side. Endoplasmic reticulum .
Tissue Specificity	Ubiquitously expressed with stronger levels detected in heart, liver and skeletal muscle.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

