



LRC19 rabbit pAb

Catalog No	YP-Ab-09305
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	LRRC19
Immunogen	Synthesized peptide derived from human LRC19 AA range: 30-80
Specificity	This antibody detects endogenous levels of human LRC19
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRC19;LRRC19;Leucine-rich repeat-containing protein 19
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in renal collecting duct epithelial cells.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

