



HDHD1 rabbit pAb

Catalog No	YP-Ab-09308
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB
Gene Name	HDHD1 DXF68S1E FAM16AX GS1 HDHD1A
Immunogen	Synthesized peptide derived from human HDHD1 AA range: 89-139
Specificity	This antibody detects endogenous levels of human HDHD1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	5'-PsiMPase;DXF68S1E;EC:3.1.3.96;FAM16AX;GS1;HDHD1;HDHD1 DXF68S1E FAM16AX GS1 HDHD1A;Pseudouridine-5'-monophosphatase (5'-PsiMPase) (EC 3.1.3.n6) (Haloacid dehalogenase-like hydrolase domain-containing protein 1) (Haloacid dehalogenase-like hydrolase domain-containing protein 1A) (Pro
Observed Band	25 kDa
Calculated Molecular Weight	228 aa, 25 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

