



CRFM7 rabbit pAb

Catalog No	YP-Ab-09311
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CHRFAM7A
Immunogen	Synthesized peptide derived from human CRFM7 AA range: 124-174
Specificity	This antibody detects endogenous levels of human CRFM7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CRFM7;CHRFAM7A;CHRNA7-FAM7A fusion protein (CHRNA7-DR1) (D-10)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in hippocampus.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

