



ST1A3 rabbit pAb

Catalog No	YP-Ab-09312
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	SULT1A3 STM; SULT1A4
Immunogen	Synthesized peptide derived from human ST1A3 AA range: 227-277
Specificity	This antibody detects endogenous levels of human ST1A3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:2.8.2.1;HAST1/HAST2;OK/SW-cl.88;P PST;P PST 1;ST1A3;SULT1A3 STM; SULT1A4;Sulfotransferase 1A3/1A4 (ST1A3/ST1A4) (EC 2.8.2.1) (Aryl sulfotransferase 1A3/1A4) (Catecholamine-sulfating phenol sulfotransferase) (HAST3) (M-PST) (Monoamine-sulfating phenol sulfotransferase) (Pla
Observed Band	34 kDa
Calculated Molecular Weight	34 kDa
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

