



HLF rabbit pAb

Catalog No	YP-Ab-09313
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	HLF
Immunogen	Synthesized peptide derived from human HLF AA range: 156-206
Specificity	This antibody detects endogenous levels of human HLF
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hepatic leukemia factor;HLF
Observed Band	33kDa, 66 kDa
Calculated Molecular Weight	295 aa, 33 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Highly expressed in liver; lower levels in lung and kidney.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

