



ZN143 rabbit pAb

Catalog No	YP-Ab-09315
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB
Gene Name	ZNF143 SBF STAF
Immunogen	Synthesized peptide derived from human ZN143 AA range: 518-568
Specificity	This antibody detects endogenous levels of human ZN143
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN143;ZNF143 SBF STAF;Zinc finger protein 143 (SPH-binding factor) (Selenocysteine tRNA gene transcription-activating factor) (hStaf)
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in all tissues tested, with the strongest expression in ovary.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

